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Letters, numerals and symbols shall be at least 12 mm high, except for packagings of 30 litres or 30 kg capacity or less, when they shall be at least 6 mm in height and for packagings of 5 liters or 5 kg or less when they shall be of an appropriate size.

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Lid
BOUILLARD SAS
38 route de Louhans
F - 71370 ST GERMAIN DU PLAIN


Test: According to EN ISO 23907 (October 2012); Sharps injury protection – Requirements and test methods – Sharps containers and NF X 30-511 (April 2015); Containers for medical care waste – Additional and/or alternative characteristics and requirements for sharps containers

Test Results: The present materials comply with the European Norm EN ISO 23907 (October 2012) and French norm NF X 30-511 (April 2015).
The complete description of the presented material, the results and findings of the tests are given in the report **PBO/pbo/P-15.176-A (F)**, **PBO/pbo/P-15.176-A (E)** and **PBO/pbo/P-15.176-B (F)** of **24-Jul-2015** which is inextricable to this certificate

The Belgium Packaging Institute Ltd, Z.1. Research Park 280, B-1731 Zellik, certifies the here above mentioned statements. Any change of the presented material or production method involves a new series of conformity tests on the modified packages to be performed by the Belgium Packaging Institute.

The certificate and report are valid for an unlimited period as far as the production method remains identical and the packaging remains within the specifications of the tested samples.

Quality Manager

Filip Van Impe

Analyst

Philippe Bourdon

page 1/1

This certificate of test and its report shall not be reproduced except in full, without written approval of the B.V.I.-I.B.E.

RECOMMENDATIONS FOR THE GOOD USE OF THE PICADOR SHARPS BOXES

The sharps boxes in polypropylene as well as their treatment are costly for your structure: they are especially designed for solid sharp infectious objects which represent a moderate risk to the community.

→ please refer to your internal waste segregation scheme as well as to the British Standard BS EN ISO 23907: 2012.

The level of security can be perceptibly deteriorated in case of misuse especially if the box is not correctly mounted or if the permanent closure system is not correctly activated.

Recommendations before use:

- Read carefully the instructions for use, especially for the assembly of the components and for the permanent closure,
- Respect the mounting instructions and precautions,
- Use the packaging only once it has been correctly mounted ie when the lid is correctly assembled with the box,
- Identify the handle and locate the filling limit,
- Attach the box on a support by using the fixing tie in order to ensure the use with one hand,
- Indicate the date at which it was put into service in the 'producer' zone provided on the sticker.


Recommendations during use:

- **Respect the filling limit: as soon as it has been reached, close permanently the box,**
- Keep out of the reach of children,
- Make sure to have used the provisory closure system after each use,
- Never introduce the hand or a part of the hand to ease the introduction of waste,
- Never force for introduction of waste,
- Check the permanent closure according to the instructions for use ,
- Never reopen the box,
- The boxes can be autoclaved at 134°C during 18 minutes prior to use. **Waste to be eliminated according the regulations.**

Box reference	Useful capacity [L]	Max. gross weight [kg]	External dimensions (L x l x h) [cm]
Pi0500			
Pi1500	1,15	1,00	18,5 x 13 x 12
Pi2500	2,13	1,38	18,5 x 13 x 20
Pi3500	2,80	1,86	22,7 x 16 x 16,8
Pi5500	4,68	2,84	22,7 x 16 x 27
Pi7500	6,38	3,83	22,7 x 16 x 39
Pi9000	7,50	4,75	18,5 x 13 x 64,5

HOW TO READ THE DATE MARKERS ON THE PICADOR

The date marker on the fixing tie of the box indicates the week of manufacturing of the box during the year (out of 52 weeks).


ex.: for instance here the box has been manufactured in week 29.

→ : this arrow points the figure in tens or the number of the week (so from 0 to 5)

? : this arrow points the unit of the number of the week (so from 0 to 9)

The date marker under the lid at the level of the push buttons indicates the month and the year of manufacturing of the lid:

← : number of the month (arrow points on it) + year (on both sides of the arrow) of manufacturing

HOW TO READ THE UN NUMBER

3H2/Y1/S/12/F/BVT319110/PLAZUR

3H2 is the code describing the type of packaging according to 6.1.2 of ADR

3: Jerrican
H: Plastic material
2: removable head

Y is for packing group II et III

1 is the indication of the maximum weight in kg

S indicates that the packaging is aimed at carrying solids

12 is the year of manufacturing of the box

,F' is for France, the state authorizing the use of the mark

BVT = name of the tests laboratory followed by the number of the UN certificate

Manufacturer of the box

The marking engraved on the bottom of the box indicates that it corresponds to a model type, which passed successfully the tests prescribed in the ADR and that it is in conformity with the Chapter 6.1.3 of the ADR, which is linked to the manufacturing (and not to the use) of the packaging.